EN Instructions for use

The Accu-Chek FastClix finger pricker is a device for simple, virtually pain-free collection of capillary blood from the fingertips. The finger pricker has 11 different penetration depths so that you can adjust the penetration depth to the texture of your skin. With the release button, you can prime the finger pricker and trigger a lancet in a single step. When a lancet is triggered, it penetrates the skin. You insert a lancet drum with 6 lancets into the finger pricker. You can obtain blood using a new, sterile lancet 6 times before the lancet drum needs changing.

These instructions for use feature the following 3 symbols:

- This symbol indicates a possible risk of injury or of damage to your health.
- This symbol draws attention to actions that could result in damage to the finger pricker.
- This symbol draws your attention to important information.

The Accu-Chek FastClix finger pricker is intended for patient self-monitoring by an individual person only. It must not be used to collect blood from more than one person as it does not incorporate any features to guard against cross-infection. There is a risk of infections being transmitted if the finger pricker is used by other people, even by family members, or if healthcare professionals use this finger pricker to obtain blood from different people. This finger pricker is therefore not suitable for professional use in healthcare facilities and institutions.

- Keep the finger pricker with all its components out of reach of children under 3 years. There is a risk of suffocation if small parts (e.g. caps or similar objects) are swallowed.

Only insert Accu-Chek FastClix lancet drums into the Accu-Chek FastClix finger pricker. Other lancets or lancet drums may impair the operation of the finger pricker or damage it permanently.

Preparation

- Wash your hands with warm water and soap before you obtain blood (see 1).
- Dry your hands thoroughly. This reduces contamination of the puncture site and stimulates blood flow.

Inserting the lancet drum

- You have to insert a lancet drum to be able to obtain blood with the finger pricker.
- Remove the cap from the finger pricker (see 2).

- Insert a new lancet drum – white end first – into the finger pricker (see 2).

The lancet drum is properly inserted when it is completely inside the finger pricker, and only the projecting edge protrudes from the finger pricker at the front (see 3).

You must not insert the lancet drum into the finger pricker and simultaneously press the release button or hold the finger pricker with the release button resting on a surface such as a table top. This could release a lancet and inadvertently cause injury.

- Place the cap back on again. The cut-out section at the edge of the cap must fit around the correspondingly shaped curvature with the penetration depth indicator (silver square) (see 4).

The cap is on correctly when it is pushed back as far as it will go, and you hear and feel it click into place.

The cap can only be correctly placed up to the stop position when the lancet drum is completely inside the finger pricker. Do not use any force to place the cap on as this might damage the cap and the finger pricker. Correct the position of the lancet drum instead.

Each lancet drum contains 6 lancets. The first lancet is now ready to be used. The window showing the number of available lancets shows 6; there are now 6 lancets available (see 5).

Do not remove the lancet drum before all 6 lancets have been used. The lancet drum cannot be re-inserted once it has been removed from the finger pricker.

Setting the penetration depth

You can set 11 penetration depths (0.5 to 5.5). The silver square (penetration depth indicator) on the finger pricker indicates the current penetration depth setting. The larger the number, the greater the penetration depth.

If you have no experience with this finger pricker, we recommend setting a low penetration depth, such as penetration depth 2.

- Rotate the cap until the desired penetration depth lines up with the indicator (see 6).

Obtaining a blood drop

In general, blood can be obtained from any finger. Certain fingers can be unsuitable if, for example, a skin or fingernail infection is present.

We recommend obtaining capillary blood from the sides of the fingertips as these areas are the least sensitive to pain.

Always use a new lancet when you obtain blood. In doing this, you reduce the risk of infection, and obtaining blood remains virtually pain free.

Only use the finger pricker when the cap is attached. Without the cap, the lancet penetrates too deeply and the puncture can become painful.

- In individual cases, the lancet may protrude from the cap after obtaining blood. There is a risk of injury or infection if you or another person touches the protruding lancet. If the lancet protrudes from the cap, carefully set the finger pricker aside. Keep the finger pricker out of reach of other people and contact the Accu-Chek Customer Careline.

- Press the finger pricker firmly against the selected puncture site (see 7).
- Press the release button all the way down (see 8).

The lancet is released and it penetrates the skin.

- Massage the finger in the direction of the fingertip to encourage a blood drop to form.

The amount of blood that emerges depends on the penetration depth and the pressure used to hold the finger pricker against the skin.

If not enough blood emerges, apply more pressure on the finger pricker when you next obtain blood.

If that is not sufficient, increase the penetration depth progressively in addition.

If too much blood emerges, decrease the penetration depth.

- Test your blood glucose immediately after you obtained blood, according to the instructions for use for your blood glucose meter.

We recommend protecting the puncture site from contamination afterwards.

Loading a new lancet

- If you have already used the finger pricker to obtain blood, you will need a new lancet the next time you obtain blood.
- Push the lever at the side of the finger pricker as far as it will go in direction A and then back in direction B (see 9).

The next lancet is now ready to be used. The window next to the lever shows the new number of available lancets, which has decreased by one, for example, 5 instead of 6.

- Proceed as in the previous step each time you need to change the lancet.

Each time a new lancet is loaded, the number of available lancets shown in the window decreases by one, for example, from 4 to 3. When you have loaded the sixth and last lancet, the window shows 1.

Replacing a used lancet drum

When you have used the sixth and last lancet, you need to replace the used lancet drum with a new one.

- Remove the cap from the finger pricker (see 2).
- Remove the lancet drum from the finger pricker (see 10).

When a lancet drum has been used, a red stripe is visible on the white cylinder.

- Insert a new lancet drum (see section *Inserting the lancet drum*).

A used lancet drum has a built-in lock that stops you from using it again.

Used lancet drums can be disposed of in household waste if no other regulations apply in your country.

Storing the finger pricker

Do not store the finger pricker at extreme temperatures, for example, in a hot car. This can impair the operation of the finger pricker.

Cleaning and disinfecting the finger pricker and cap

To prevent infections from being transmitted, you must clean and disinfect the finger pricker and cap regularly:


- At least once a week
- When there is blood on the finger pricker or cap
- Always before someone else handles the finger pricker, for example, to assist you.

You need lint-free cloths in which you can wrap the whole finger pricker, water and 70 % isopropanol. Before you disinfect the finger pricker and the cap, you first have to clean them to remove blood and other dirt.

To avoid damaging the finger pricker and the cap:

- Do not dip the finger pricker and the cap in the fluids.
- Use only water for cleaning and only 70 % isopropanol for disinfecting.

- Remove the cap from the finger pricker (see 2).

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