High Viscosity
DERMABOND™

TOPICAL SKIN ADHESIVE

(2-Octyl Cyanoacrylate)

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ETHICON, LLC
San Lorenzo, Puerto Rico 00754
© ETHICON, Inc. 2009

High Viscosity **DERMABOND™**

Topical Skin Adhesive (2-Octyl Cyanoacrylate)

DESCRIPTION

DERMABOND™ Topical Skin Adhesive is a sterile, liquid topical skin adhesive containing a monomeric (2-octyl cyanoacrylate) formulation and the colorant D & C Violet #2. It is provided in a single use applicator packaged in a blister pouch. The applicator is comprised of a crushable glass ampule contained within a plastic vial with attached applicator tip. The product is also available with a pen applicator. As applied to the skin, the liquid is syrup-like in viscosity and polymerizes within minutes. Studies have shown that following application of DERMABOND, it acts as a barrier to prevent microbial infiltration into the healing wound. See DIRECTIONS FOR USE.

INDICATIONS

DERMABOND Topical Skin Adhesive is intended for topical application only to hold closed easily approximated skin edges of wounds from surgical incisions, including punctures from minimally invasive surgery, and simple, thoroughly cleansed, trauma-induced lacerations. DERMABOND adhesive may be used in conjunction with, but not in place of, deep dermal stitches.

CONTRAINDICATIONS

- Do not use on any wound with evidence of active infection, gangrene, or wounds of decubitus etiology.
- Do not use on mucosal surfaces or across mucocutaneous junctions (e.g., oral cavity, lips), or on skin which may be regularly exposed to body fluids or with dense natural hair.
- Do not use on patients with a known hypersensitivity to cyanoacrylate or formaldehyde.

WARNINGS

- DERMABOND adhesive is a fast setting adhesive capable of adhering to most body tissue and many other materials, such as latex gloves and stainless steel. Inadvertent contact with any body tissue, and any surfaces or equipment that are not disposable or that cannot be readily cleaned with a solvent such as acetone should be avoided.
- Polymerization of DERMABOND adhesive may be accelerated by water or fluids containing alcohol: DERMABOND adhesive should not be applied to wet wounds.
- DERMABOND adhesive should not be applied to the eye. If contact with the eye occurs, flush the eye copiously with saline or water. If residual adhesive remains, apply topical ophthalmic ointment to help loosen the bond and contact an ophthalmologist.
- When closing facial wounds near the eye with DERMABOND adhesive, position the patient so that any run-off of adhesive is away from the eye. The eye should be closed and protected with gauze. Prophylactic placement of petroleum jelly around the eye, to act as a mechanical barrier or dam, can be effective in preventing inadvertent flow of adhesive into the eye. DERMABOND adhesive will not adhere to skin pre-coated with petroleum jelly. Therefore, avoid using petroleum jelly on any skin area where DERMABOND adhesive is intended to adhere. Use of DERMABOND adhesive near the eye has inadvertently caused some patient's eyelids to be sealed shut. In some of these cases, general anesthesia and surgical removal has been required to open the eyelid.
- DERMABOND adhesive should not be used below the skin because the polymerized material is not absorbed by tissue and can elicit a foreign body reaction.
- DERMABOND adhesive should not be used in high skin tension areas or across areas of increased skin tension, such as knuckles, elbows, or knees, unless the joint will be immobilized during the skin healing period or unless skin tension has been removed by application of another wound closure device (e.g., sutures or skin staples) prior to application of high viscosity DERMABOND adhesive.
- DERMABOND adhesive treated wounds should be monitored for signs of infection. Wounds with signs of infection, such as erythema, edema, warmth, pain and pus, should be evaluated and treated according to standard practice for infection.

- DERMABOND adhesive should not be used on wound sites that will be subjected to repeated or prolonged moisture or friction.
- DERMABOND adhesive should only be used after wounds have been cleaned and debrided in accordance with standard surgical practice. Local anesthetic should be used when necessary to assure adequate cleansing and debridement.
- Excessive pressure of the applicator tip against wound edges or surrounding skin can force the wound edges apart and allow adhesive into the wound. Adhesive within the wound could delay wound healing and/or result in adverse cosmetic outcome. Therefore, DERMABOND adhesive should be applied with a very light brushing motion of the applicator tip over easily approximated wound edges.
- DERMABOND adhesive polymerizes through an exothermic reaction in which a small amount of heat is released. With the proper technique of applying DERMABOND adhesive in multiple thin layers (at least 2) onto a dry wound and allowing time for polymerization between applications, heat is released slowly and the sensation of heat or pain experienced by the patient is minimized. However, if DERMABOND adhesive is applied so that large droplets of liquid are allowed to remain unspread, the patient may experience a sensation of heat or discomfort.
- DERMABOND adhesive is packaged for single patient use. Discard remaining opened material after each wound closure procedure.
- Do not resterilize DERMABOND adhesive.
- Do not place DERMABOND adhesive in a procedure pack/tray that is to be sterilized prior to use. Exposure of DERMABOND adhesive, after its final manufacture, to excessive heat (as in autoclaves or ethylene oxide sterilization) or radiation (such as gamma or electron beam), is known to increase its viscosity and may render the product unusable.

PRECAUTIONS

- Do not apply liquid or ointment medications or other substances to the wound after closure with DERMABOND adhesive, as these substances can weaken the polymerized film and allow for wound dehiscence. DERMABOND adhesive permeability by topical medications has not been studied.
- DERMABOND adhesive permeability by fluids is not known and has not been studied.
- DERMABOND adhesive, as a liquid, is syrup-like in viscosity. To prevent inadvertent flow of liquid DERMABOND adhesive to unintended areas: (1) the wound should be held in a horizontal position, with DERMABOND adhesive applied from above, and (2) DERMABOND adhesive should be applied in multiple (at least 2), thin layers rather than in a few large droplets.
- Hold applicator away from yourself and the patient and break ampule close to its center one time only. Do not crush the contents of the applicator tube repeatedly as further manipulation of the applicator may cause glass shard penetration of the outer tube.
- DERMABOND adhesive should be used immediately after crushing the glass ampule as the liquid adhesive will not flow freely from the applicator tip after a few minutes.
- If unintended bonding of intact skin occurs, peel, but do not pull the skin apart. Petroleum jelly or acetone may help loosen the bond. Other agents such as water, saline, Betadine® Antibiotics, HIBICLENS® (chlorhexidine gluconate), or soap, are not expected to immediately loosen the bond.
- Safety and effectiveness of DERMABOND adhesive on wounds of patients with peripheral vascular disease, insulin dependent diabetes mellitus, blood clotting disorders, personal or family history of keloid formation or hypertrophy, or burst stellate lacerations, have not been studied.
- Safety and effectiveness of DERMABOND adhesive on the following wounds have not been studied: animal or human bites, puncture or stab wounds.
- Safety and effectiveness on wounds that have been treated with DERMABOND adhesive and then exposed for prolonged periods to direct sunlight or tanning lamps have not been studied.
- Safety and effectiveness of DERMABOND adhesive on wounds in vermilion surfaces has not been studied.

ADVERSE REACTIONS

- Adverse reactions encountered during clinical study:

<i>Clinical Study Outcomes</i>	<i>No Subcuticular Sutures</i>		<i>With Subcuticular Sutures</i>	
	DERMABOND	Control	DERMABOND	Control
	N (%)	N (%)	N (%)	N (%)
Accounting				
N, patients enrolled	240	243	167	168
N, patients treated	239	242	167	166
Patients completed	228 (95%)	215 (88%)	164 (98%)	162 (96%)
Adverse Reactions				
Suspected Infection*	8 (3.6%)	2 (0.9%)	6 (3.6%)	2 (1.2%)
Wound type				
# Lacerations	8	2	1	0
# Incisions	0	0	5	2
Dehiscence with Need for Retreatment	6 (2.5%)	5 (2.1%)	3 (1.8%)	0
Acute Inflammation				
Erythema	26 (11.5%)	74 (33.0%)	52 (31.3%)	75 (45.1%)
Edema	22 (9.7%)	28 (12.5%)	62 (37.3%)	71 (42.8%)
Pain	14 (6.1%)	13 (5.8%)	56 (33.7%)	57 (34.3%)
Warmth	3 (1.3%)	6 (2.6%)	3 (1.8%)	4 (2.4%)

*In the clinical study, presence of infection was to be identified by observation of redness more than 3–5 mm from the repaired wound, swelling, purulent discharge, pain, increased skin temperature, fever, or other systemic signs of infection. (See clinical study). Confirmatory culture was not routinely obtained. Among cases of suspected infection for DERMABOND adhesive, 7/14 (50%) were in patients less than 12 years old with traumatic lacerations; overall, 8 of the 14 (approximately 60%) DERMABOND adhesive wounds with suspected infections were associated with sub-optimal cosmetic outcome.

- Reactions may occur in patients who are hypersensitive to cyanoacrylate or formaldehyde. See CONTRAINDICATIONS.
- The polymerization of DERMABOND adhesive on the skin releases small amounts of heat which may cause a sensation of heat or discomfort in some patients.
- Adverse reactions may be experienced following DERMABOND adhesive contact with the eye.

CLINICAL STUDIES

Clinical Study Comparing High Viscosity DERMABOND Adhesive and Low Viscosity DERMABOND Adhesive for Closure of Trauma-Induced Lacerations

Description: A prospective, randomized, controlled, unmasked study was conducted to evaluate the safety and effectiveness of closing the approximated skin edges of trauma-induced lacerations using high viscosity DERMABOND adhesive in comparison to the currently marketed low viscosity DERMABOND adhesive, with or without stitches placed below the skin surface according to investigator judgment.

The study population included patients at least one year of age, in good general health, who signed informed consent and agreed to follow-up visits. Patients were excluded if presenting with: significant multiple trauma, peripheral vascular disease, insulin dependent diabetes mellitus, blood clotting disorder, keloid formation or hypertrophy history (patient or family), cyanoacrylate or formaldehyde allergy, burst or stellate lacerations due to crush or hard blow, animal or human bite, and decubitus ulcer. One unit of either the high viscosity or low viscosity DERMABOND product was used to close the wound. Wound length and width were measured in millimeters; wound depth was not measured. Most wounds in the study were clean, small, superficial lacerations which did not penetrate the dermis completely nor sufficiently to require dermal suture placement (average wound length=18.7mm).

Follow-up was at 14 days (± 2 days) and at 30 days (± 2 days). The Modified Hollander Cosmesis Scale (MHCS), a validated scale, was used to evaluate cosmesis at the 30-day (± 2 days) follow-up visit. This scale evaluates step-off borders, edge inversion, contour irregularities, excess inflammation, wound margin separation, and overall wound appearance.

Results: The primary measure of device effectiveness in the study was wound closure at day 14, defined as continuous approximation of wound margins from the time of wound closure until the day of evaluation. Results indicate that high viscosity DERMABOND adhesive was equivalent to the low viscosity DERMABOND adhesive control for effectiveness of wound closure at day 14. Secondary effectiveness measures included an assessment of migration of the liquid adhesive from the wound site during application and an assessment of the presence of the polymer film on the wound at the time of the 14 day follow-up. Results show a significant reduction in the occurrence of migration of the liquid adhesive from the wound site during application for the high viscosity DERMABOND adhesive compared to the low viscosity DERMABOND adhesive control. No significant difference was observed between the two treatment groups for the presence of the polymer film at day 14.

Clinical Study Comparing Low Viscosity DERMABOND Adhesive to Sutures, Staples, and Adhesive Strips

Description: A prospective, randomized, controlled, unmasked study was conducted to evaluate the safety and effectiveness of closing the approximated skin edges of surgical incisions, including punctures from minimally invasive surgery, and trauma-induced lacerations using low viscosity DERMABOND adhesive in comparison to U.S.P. size 5-0 or smaller sutures, adhesive strips or staples, with or without dermal closure (subcuticular stitch) as per investigator judgment.

**Summary of Effectiveness Results Comparing DERMABOND Adhesive to Sutures
(U.S.P. size 5-0 and smaller diameter), Staples, and Adhesive Strips**

<i>Clinical Study Outcomes</i>	<i>NSS</i>		<i>WSS</i>	
	DERMABOND	Control	DERMABOND	Control
	N (%)	N (%)	N (%)	N (%)
<i>Accounting</i>				
N, patients enrolled	240	243	167	168
N, patients treated	239	242	167	166
Patients completed	228 (95%)	215 (88%)	164 (98%)	162 (96%)
N, control: suture/strips/staples/missing		194/46/1/1		116/45/5/0
<i>Wound Closure Assessment</i>				
Immediate: Additional Devices	18 (7.5%)	13 (5.4%)	2 (1.2%)	11 (6.6%)
@ 5–10 days: 100% epidermal apposition	169 (75.1%)	199 (88.8%)	140 (84.3%)	160 (96.4%)
>50% epidermal apposition	205 (91.1%)	214 (95.5%)	163 (98.2%)	165 (99.4%)
@ 3 months: Cosmesis Score* = 0 (optimal)	188 (82.5%)	180 (83.7%)	128 (78.0%)	128 (79.0%)
Median Time for Treatment (Minutes)	1.5	6.0	1.3	2.9

*** Cosmesis: Modified Hollander Cosmesis Scale**

The study population included patients at least one year of age, in good general health, who signed informed consent and agreed to follow-up visits. Patients were excluded if presenting with: significant multiple trauma, peripheral vascular disease, insulin dependent diabetes mellitus, blood clotting disorder, keloid formation or hypertrophy history (patient or family), cyanoacrylate or formaldehyde allergy, burst or stellate lacerations due to crush or hard blow, animal or human bite, and decubitus ulcer.

Follow-up was at 5–10 days and at 3 months. All wounds were assessed by visual inspection at 5–10 days after wound closure. The total kinds of wounds treated in the study were 46.1% lacerations and 53.9% incisions. The incisions were comprised of 47.8% excisions of skin lesions, 27.4% minimally invasive surgery punctures, and 24.8% general surgery incisions.

For wounds closed without subcuticular stitches, mean wound length was 1.5cm, mean wound width was 2.5mm, and mean wound depth was 5.8mm. For wounds closed with subcuticular stitches, mean wound length was 3.2cm, mean wound width was 5.3mm, and mean wound depth was 3.8mm.

If the primary method of closure was insufficient for closure, an additional securing device was placed. The time to perform treatment included the time required later to remove the closure device when applicable.

The Modified Hollander Cosmesis Scale (MHCS), a validated scale, was used to evaluate cosmesis at three months: step-off borders, edge inversion, contour irregularities, excess inflammation, wound margin separation, and overall appearance.

DIRECTIONS FOR USE

1. The application of DERMABOND adhesive requires thorough wound cleansing. Follow standard surgical practice for wound preparation before application of DERMABOND adhesive (i.e., anesthetize, irrigate, debride, obtain hemostasis and close deep layers).
2. Pat the wound dry with dry, sterile gauze to assure direct tissue contact for adherence of the DERMABOND adhesive to the skin. Moisture accelerates DERMABOND adhesive's polymerization and may affect wound closure results.
3. To prevent inadvertent flow of liquid DERMABOND adhesive to unintended areas of the body, the wound should be held in a horizontal position and the DERMABOND adhesive should be applied from above the wound.
4. DERMABOND adhesive should be used immediately after crushing the glass ampule, since the liquid DERMABOND adhesive will flow freely from the tip for only a few minutes. Remove the applicator from the blister pouch. If using the pen applicator, refer to the instructions on the pouch for crushing the glass ampule and expressing the liquid adhesive. If using the plastic vial, hold the applicator with the thumb and a finger and away from the patient to prevent any unintentional placement of the liquid DERMABOND adhesive into the wound or on the patient. While holding the applicator, and with applicator tip pointed upward, apply pressure at the midpoint of the ampule to crush the inner glass ampule. Invert and gently squeeze the applicator just sufficiently to express the liquid DERMABOND adhesive to moisten the applicator tip.
5. Approximate wound edges with gloved fingers or sterile forceps. Slowly apply the liquid DERMABOND adhesive in multiple (at least 2) thin layers to the surface of the approximated wound edges using a gentle brushing motion. Wait approximately 30 seconds between applications or layers. Maintain manual approximation of the wound edges for approximately 60 seconds after the final layer.

NOTE: DERMABOND adhesive polymerizes through an exothermic reaction. If the liquid DERMABOND adhesive is applied so that large droplets are allowed to remain without being evenly spread, the patient may experience a sensation of heat or discomfort. The sensation may be higher on sensitive tissues. This can be minimized by applying DERMABOND adhesive in multiple thin layers (at least 2).

NOTE: Excessive pressure of the applicator tip against the wound edges or surrounding skin can result in forcing the wound edges apart and allowing DERMABOND adhesive into the wound. DERMABOND adhesive within the wound could delay wound healing and/or result in adverse cosmetic outcome.

NOTE: Full apposition strength is expected to be achieved about 2.5 minutes after the final layer is applied, although the top adhesive layer may remain tacky for up to approximately 5 minutes. Full polymerization is expected when the top DERMABOND adhesive layer is no longer sticky.

6. Do not apply liquid or ointment medications onto wounds closed with DERMABOND adhesive because these substances can weaken the polymerized film, leading to wound dehiscence.
7. Protective dry dressings such as gauze, may be applied only after DERMABOND adhesive film is completely solid/polymerized: not tacky to the touch (approximately five minutes after application). Allow the top layer to fully polymerize before applying a bandage.
If a dressing, bandage, adhesive backing or tape is applied before complete polymerization, the dressing can adhere to the film. The film can be disrupted from the skin when the dressing is removed, and wound dehiscence can occur.
8. Patients should be instructed to not pick at the polymerized film of DERMABOND adhesive. Picking at the film can disrupt its adhesion to the skin and cause dehiscence of the wound. Picking at the film can be discouraged by an overlying dressing.
9. Apply a dry protective dressing for children or other patients who may not be able to follow instructions for proper wound care.
10. Patients treated with DERMABOND adhesive should be provided the printed instruction sheet entitled, How to Care for Your Wound After It's Treated With DERMABOND Topical Skin Adhesive. This instruction sheet should be reviewed with each patient or guardian to assure understanding of the proper care for the treatment site.

- Patients should be instructed that until the polymerized film of DERMABOND adhesive has sloughed naturally (usually in 5–10 days), there should be only transient wetting of the treatment site. Patients may shower and bathe the site gently. The site should not be scrubbed, soaked, or exposed to prolonged wetness until after the film has sloughed naturally and the wound has healed closed. Patients should be instructed not to go swimming during this period.
- If removal of DERMABOND adhesive is necessary for any reason, carefully apply petroleum jelly or acetone to the DERMABOND film to help loosen the bond. Peel off the film, do not pull the skin apart.

HOW SUPPLIED

DERMABOND adhesive is supplied sterile, in a pre-filled, single-use applicator. The applicator is comprised of a crushable glass ampule contained within a plastic vial with attached applicator tip. This product is also available with a pen applicator. The applicator contains the liquid adhesive. The applicator is packaged in a blister pouch to maintain the device sterile until opened or damaged.

DERMABOND adhesive may be available in boxes of 6 or 12 applicators.

STORAGE

Recommended storage conditions: below 30°C, 86°F, away from moisture, direct heat, and direct light. Do not use after expiry date.

STERILITY

DERMABOND adhesive is originally sterilized by dry heat and ethylene oxide gas. Do not resterilize. Do not use if package is opened or damaged. Discard any unused material following completion of medical procedure.

STERILE SINGLE USE ONLY

REPORTING

Physicians should use the following toll free number 1-800-255-2500 (valid in U.S.A. only), when reporting adverse reactions or potentially threatening complications involving DERMABOND adhesive.

CAUTION

Federal (U.S.A.) Law restricts this device to sale by or on the order of a physician.

©Trademark of Purdue Fredericks

†Registered Trademark of Mölnlycke RM Ltd.

 <p>0086 CE mark and identification number of Notified Body. Product conforms to the essential requirements of the Medical Device Directive 93/42/EEC</p>	 <p>Light Sensitive</p>	 <p>30°C 86°F max</p>	
 <p>Use by — year & month</p>	<p>STERILE ↓</p> <p>Method of Sterilization — Dry Heat</p>	<p>Store at room temperature. (Below 30°C, 86°F)</p>	<p>Contents: Vial</p>
<p>LOT</p> <p>Batch number</p>	<p>STERILE EO</p> <p>Method of Sterilization — Ethylene Oxide</p>	 <p>See instructions for use</p>	
	 <p>Do not resterilize</p>	 <p>Do not reuse</p>	<p>Crush Vial After Opening the Package</p>

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Johnson & Johnson Medical Limited
ETHICON S.A.S.
Medical Limited
JOHNSON & JOHNSON MEDICAL S.p.A.
CH: Johnson & Johnson AG