ACCU-CHEK and ACCU-CHEK FLEXLINK
are trademarks of Roche.

© 2016 Roche Diabetes Care

Roche Diabetes Care GmbH
Sandhofer Strasse 116
68305 Mannheim, Germany
www.accu-chek.com

Last update: 2016-03

04537246001/H

ACCU-CHEK[®] FlexLink

Infusion set

PDF

Approved according to
SOP 04.07.08 – Anl. A, Vers. 02

Kat. 03


The Signature will not be printed!

LAY


LAB

LAN


1


2


3


4


5


6

7

8


Instructions for use	4	Návod k použití	88
Gebrauchsanweisung	16	Instrukcja obsługi	100
Instructions d'utilisation	28	Brugsanvisning	112
Istruzioni per l'uso	40	Bruksanvisning	124
Instrucciones de uso	52	Bruksanvisning	136
Gebruiksaanwijzing	64	Käyttöohje	148
Instruções de utilização	76		

Instructions for use	en
Gebrauchsanweisung	de
Instructions d'utilisation	fr
Istruzioni per l'uso	it
Instrucciones de uso	es
Gebruiksaanwijzing	nl
Instruções de utilização	pt
Návod k použití	cs
Instrukcja obsługi	pl
Brugsanvisning	da
Bruksanvisning	sv
Bruksanvisning	no
Käyttöohje	fi

ACCU-CHEK[®]
FlexLink

Infusion set

Indications for use

The Accu-Chek FlexLink infusion set is intended for the subcutaneous infusion of insulin administered with microdosage insulin pumps.

Contents of the infusion set packages

Package 1: head set and transfer set

- A** introducer needle handle
- B** protective cover
- C** introducer needle
- D** soft cannula
- E** adhesive backing
- F** self-adhesive
- G** cannula housing

- H** connector
- I** side clips of the connector
- J** tubing
- K** luer-lock connector
- L** disconnect cover
- M** blue holder cap

Package 2: head set

Components A–G, L and M from package 1.
May only be used together with the transfer set from package 1.

en

⚠ WARNINGS

Read the instructions carefully before use. We recommend that you consult your insulin pump manual for information regarding pump therapy.

- When using the infusion set for the first time, your doctor or healthcare team must be present.
- This device remains sterile until the package is opened or becomes damaged. Do not use if the sterile package has been opened or damaged.
- Protect the product from excessive humidity, sunlight and heat. Store at room temperature.
- Dispose of your infusion set after use

according to local regulations. Do not clean or re-sterilize.

- The infusion set is for subcutaneous use only.
- Always screw the luer-lock connector firmly into the pump adapter. If this is not done, a loose connection between the infusion set and the pump could cause leakage. Avoid using tools to screw the luer-lock connector into place, otherwise damage could result to the luer-lock connector.
- Do not leave air bubbles in the infusion set. Fill carefully.
- A small percentage of Accu-Chek FlexLink infusion sets using soft cannulas can accidentally become crimped during insertion or displaced while in use. If

this should occur, use a new infusion set immediately.

- Check the infusion site frequently to make sure the soft cannula remains firmly in place. Because the cannula is so soft, you will not notice any pain if it is pulled out. The soft cannula must always be completely inserted in order to receive the full amount of insulin. Replace the infusion set as soon as the self-adhesive loosens.
- Do not reinsert the introducer needle into the soft cannula. Reinsertion could cause tearing of the soft cannula which would result in unpredictable insulin delivery.
- If you have to disconnect your infusion set pay particular attention to hygiene.
- Never use a soft cannula for more than

72 hours (3 days). Never use the tubing for more than 6 days. When changing the cartridge, make sure that the transfer set is completely filled before connecting it to the head set.

- If the infusion site becomes inflamed, replace your infusion set and use a new site.
- Do not put disinfectants, perfumes, deodorants or insect repellent on your infusion set as these may affect the integrity of the infusion set.
- Never prime the transfer set or attempt to free a clogged line while the transfer set is connected to the head set. You may accidentally infuse an uncontrolled quantity of insulin.
- Avoid mechanical stress on the infusion site

en

and on the infusion set. Pay attention when carrying heavy weights.

- Before inserting the infusion set, prepare the infusion site in accordance with the instructions of your healthcare professional.
- Check your blood glucose level 1 to 3 hours after inserting your infusion set and inspect the infusion site on a regular basis. Do not insert a new soft cannula just prior to bedtime. Check your blood glucose frequently (check with your healthcare professional).
- Check your infusion set at the luer-lock connector to tube connection periodically (every 3 hours during the day and before going to sleep) during use. A break may result in leakage, causing interruption of

insulin delivery. If you notice a break, change your infusion set immediately and check your blood glucose.

- If there is an unexplained increase in your blood glucose level, causing it to be too high, or if an “occlusion” alarm occurs, check for clogs and leaks. If in doubt, change your infusion set.
- This product contains small parts. There is a risk of suffocation if small parts (e.g. covers, caps or similar objects) are swallowed.
- Pay attention to the tubing of the infusion set connected to your body. There is a risk of strangulation if the tubing becomes wrapped around your neck. In case of doubt, use the shortest tubing length to minimise the risk of strangulation.

Insert the infusion set

Picture 1

Choose a site away from the waistline, bones, scar tissue, navel and recent infusion sites. Wash your hands and clean the infusion site according to the instructions from your healthcare professional. Allow the infusion site to dry before the soft cannula is inserted.

Picture 2

While holding the infusion set at the blue holder cap (M) with thumb and forefinger, remove the transparent protective cover (B) by turning and pulling. Keep the protective cover for later use in order to safely dispose of the introducer needle (C).

Picture 3

Remove both pieces of the adhesive backing (E) from the self-adhesive (F).

Picture 4

Pinch the subcutaneous tissue. Insert the introducer needle (C) and soft cannula (D) quickly and straight into the infusion site, at a 90° angle. Watch the soft cannula while inserting to assure that it goes in smoothly without kinking or crimping.

en

⚠ WARNING

Ensure clean application. Do not touch the soft cannula or introducer needle. Make sure that the soft cannula is fully inserted into the subcutaneous tissue and not kinked at the skin surface.

Once the introducer needle punctures the skin, the soft cannula must be quickly pushed in the rest of the way, or it may kink upon insertion, causing reduced or blocked insulin delivery.

Picture 5

Carefully smooth the self-adhesive (F) onto the skin and run your fingers over the self-adhesive to ensure good contact with the skin.

Picture 6

Hold the self-adhesive (F) in place, gently press the side clips of the blue holder cap (M) and pull it from the cannula housing (G).

Picture 7

To remove the introducer needle (C) hold the cannula housing (G) firmly to your body. Place one finger on the top of the cannula housing and withdraw the introducer needle in a straight manner using the introducer needle handle (A).

Remember to dispose of the introducer needle (C) safely (needle container). Dispose of your introducer needle according to local regulations. If you inserted the infusion set by yourself and no safe needle container is at hand, consider using the protective cover to cover the needle. Ensure

that no one can injure themselves on the needle to avoid risk of infection.

Picture 8

Attach the luer-lock connector (K) to the pump by screwing it firmly into the pump adapter. Follow the pump manufacturer's User's Manual to prime the transfer set. Prime the transfer set until insulin drips from the connector needle tip. Do not leave air bubbles in either the cartridge or the transfer set. Connect the transfer set to the newly inserted head set and fill the empty space of the cannula housing with a bolus of 1.0 units (U100 insulin).

⚠ WARNING

You must fill the empty space of the newly inserted head set after connecting the transfer set. Failure to do so will result in missed insulin.

If the coupling does not click noticeably and audibly together, use a new head set and a new transfer set.

Filling volume of transfer set:

12"/ 30 cm: approx. 6 units* approx. 60 µl
24"/ 60 cm: approx. 11 units* approx. 110 µl
31"/ 80 cm: approx. 14 units* approx. 140 µl
43"/ 110 cm: approx. 18 units* approx. 180 µl

en

Filling volume of head set:

approx. 1.0 units* approx. 10 μ l

* U100 Insulin

 WARNING

Change the head set every 2 to 3 days (max. 72 hours), otherwise you risk insulin resistance or inflammation of the infusion site. Never use the transfer set for more than 6 days.

Disconnect the infusion set

en

WARNING

Consult your healthcare professional on how to compensate for any missed insulin while you are disconnected. Carefully monitor your blood glucose level while disconnected from the pump and after you reconnect.

The infusion set allows you to temporarily disconnect from your pump for showering, swimming, etc., without changing it.

Hold the self-adhesive (F) in place, gently press the side clips of the connector (I) and pull the connector (H) from the cannula housing (G). Put the disconnect cover (L) on the cannula housing.

Reconnect the infusion set

WARNING

Before reconnecting, ensure that a drop of insulin appears at the tip of the connector needle.

Fill the infusion set until insulin comes out of the connector needle. Hold the self-adhesive (F) in place and remove the disconnect cover (L) from the cannula housing (G). Re-attach the connector (H) to the cannula housing. An audible “click” ensures the connector is locked.

WARNING

It is not necessary to give an additional bolus.

Australia

Accu-Chek Enquiry Line: 1800 251 816

Pump Support: 1800 633 457

www.accu-chek.com.au

United Kingdom

Distributed in the United Kingdom by:

Roche Diabetes Care Limited

Charles Avenue, Burgess Hill

West Sussex, RH15 9RY, **United Kingdom**

Accu-Chek Pump Careline ¹⁾:

UK Freephone number: 0800 731 22 91

ROI Freephone number: 1 800 88 23 51

¹⁾ calls may be recorded for training purposes
Some mobile operators may charge for calls to these numbers.

burgesshill.insulinpumps@roche.com

www.accu-chek.co.uk

www.accu-chek.ie

New Zealand

Accu-Chek Enquiry Line: 0800 80 22 99

Pump Hotline: 0800 696 696

www.accu-chek.co.nz


Roche Diabetes Care South Africa (Pty) Ltd

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