INSTRUCTIONS FOR USE

PLEASE REFER TO THE INSTRUCTIONS FOR USE

THESE INSTRUCTIONS ARE IMPORTANT. PLEASE KEEP THEM FOR FUTURE REFERENCE.

Dear customer,

Thank you for choosing Pic Solutions easY RAPID, an automatic blood pressure monitor designed and manufactured according the most up-to-date technologies. Before using the device, we highly recommend that you read the short list of warnings contained in this manual, to make sure you have correctly understood how to use it.

INTRODUCTION

The regular use of a Blood Pressure Monitor allows you to monitor blood pressure continuously, but it should not take the place of regular check-ups by your doctor. We recommend that you visit your doctor regularly for general check-ups and for more detailed information on blood pressure. Blood pressure measurements determined with this unit are equivalent to those obtained by a trained observer using cuff/ stethoscope auscultation method, within the limits prescribed by the international ESH Validation protocol.


WARRANTY

The product is guaranteed for 5 years against all conformity defects under the normal utilisation conditions provided in the user instructions. Therefore, the warranty shall not apply in case of damages caused by improper use, wear, or accidental events. Batteries, cuff and wearing parts are not included in the warranty. In case of malfunction, contact exclusively the authorised Pic Technical Customer Service.

RAPID-TECH Comfort Experience™

This unit uses the oscillometric method to detect your blood pressure. RAPID-TECH Comfort Experience™ is the latest generation technology chosen by Pic for the blood pressure monitors of the line RAPID. It allows fast and, therefore, less annoying monitoring because the measurement is performed already during the cuff inflation phase. Pic easY RAPID has been clinically validated in accordance with the ESH (European Society of Hypertension) protocol, as a guarantee of its accuracy and precision.

NAME/FUNCTION OF EACH PART

List of included items: Blood pressure monitor , M-L size cuff (REF 0000782700000), 4xAAA batteries, user manuals, carrying bag. Features: 60x41 mm LCD display, 60 recorded measurements, Inflation-type measurement.

ABOUT BLOOD PRESSURE

What is blood pressure?

Blood pressure is the pressure exerted on the artery tube while blood flows through the arteries. The pressure measured when the heart contracts and sends blood out of the heart is systolic (highest) blood pressure. The pressure measured when the heart dilates with blood flowing back into the heart is called diastolic (lowest) blood pressure.

Why measure your blood pressure?


Among the various health problems afflicting modern people, problems associated with high blood pressure are by far the most common. High blood pressure's dangerously strong correlation with cardiovascular diseases and high morbidity has made measuring blood pressure a necessity of identifying those at risk.

Blood Pressure Fluctuation

Blood pressure fluctuates all the time! You should not be over worried if you encountered two or three measurements at high levels. Blood pressure changes over the month and even throughout the day. It is also influenced by the season, temperature and activities.

Classification of blood pressure

The blood pressure classification published by World Health Organization (WHO) and International Society of Hypertension (ISH) in 1999 is as follows:

Note: It is not intended to provide a basis of any type of rush toward emergency conditions/diagnosis based on the scheme: the scheme is meant only to discriminate between the different levels of blood pressure.

WARNING
Only a physician can tell your normal BP range. Please contact a physician if your measuring result falls out of the range. Kindly note that only a physician could tell whether your blood pressure value has reached a dangerous point.

Irrregular Heartbeat Detector

This blood pressure monitor is equipped with an algorithm for Irregular Heartbeat (IHB) Detection. During each measurement, this algorithm records the heartbeat intervals and works out the standard deviation. If the calculated value is larger than or equal to 15, this algorithm will light up the IHB symbol on the screen when displaying the measuring result.

WARNING
The appearance of the IHB icon indicates that a pulse irregularity consistent with an irregular heartbeat was detected during measurement. Usually this is NOT a cause for concern. However, if the symbol appears often, we recommend you seek medical advice. We wish to underline that the device does not replace a visit to the cardiologist, but it is useful to detect any irregularities in advance.

Why the blood pressure I get from the hospital is different from home?


The blood pressure is different even during 24 hour because of the weather, emotion, exercise etc, specially the "white coat" in hospital which makes the results are higher than the ones at home.

If the result is the same if measuring on the right arm?

Measurements can be taken on either arm; results can however differ sometimes so it is always advisable to use the same arm.

LCD display explanation

At starting check the completeness of the LCD segments

SYMBOL	DESCRIPTION	EXPLANATION
SYS	Systolic pressure	Maximum/systolic pressure result
DIA	Diastolic pressure	Minimum/diastolic pressure result
Pul/min	Pulse	Pulse/minute
▼	Deflating	The arm cuff is deflating
88	Memory	If "M" shows, the displayed measurement values is from the memory
mmHg	mmHg	Blood pressure unit of measurement mmHg
Lo +	Low battery	Batteries are low and need to be replaced
Heartbeat	Irregular heartbeat	Irregular heartbeat
Grade	Blood pressure level indicator	
MM/DD	Date	"MM" shows the month, "DD" shows the day
Avg	Average	The average of blood pressure

The Choice of Power Supply

1. Battery powered mode: 4 batteries 1,5V AAA; 6V DC.
2. Power supply via mains adapter: Input: 100-240V~, 50-60Hz, 400mA - Output: 6V 1A
- Can be powered only with the "Adapter for blood pressure monitor" REF 0000783000000 (KH0601000EW-PP) with European plug and the "Adapter for blood pressure monitor UK" REF 0000764800000 (KH0601000BW-PP) with UK plug (Not included).
- Please unplug the adaptor to depart from the using utility power.

NB: The adaptor interface is located on the right side of the monitor. Place NO obstacles on the right side for easy pull-off adaptor.

WARNING

In order to protect your device and yourself and to get the best performances, use only the batteries and adaptors according to the information provided in this manual (4 batteries 1,5V AAA; 6V DC).

Installing and Replacing the Batteries

1. Slide off the battery cover.
2. Install the batteries by matching the correct polarity, as shown.
3. Replace the cover.

Replace the batteries whenever the below happen

Is displayed **Lo +**, the display dims, the display does not light up

WARNING

- If the device is not used for a while, please remove the batteries.
- The old battery is harmful to the environment, so please disposal accordingly to the indication reported in this user manual (EU Directive 2006/66 CE).

- Do not dispose of batteries in fire. Batteries may explode or leak.
- For instructions on battery disposal please refer to the indication reported in this user manual.
- There are no use serviceable parts inside. Batteries or damage from old batteries is not covered by warranty.
- Only use brand batteries. Always replace with new batteries together. Use batteries of the same brand and same type. Other types of batteries can cause explosions or fires.
- The battery can explode if handled carelessly. Never recharge, remove or position the batteries near sources of heat. Do not incinerate.
- Rechargeable batteries are not recommended.

Setting Date and Time

It is important to set the clock before using your blood pressure monitor, so that a time stamp can be assigned to each measurement that is stored in the memory. (year: 2015–2050, time: 12H/24H). In the setting mode, use the "MEM" button to change the value and the "SET" button to confirm your choice.

1. When the unit is off, hold pressing "SET" for 3 seconds to enter the mode for year setting.

2. Press the "MEM" to change the [YEAR].

3. When you get the right year, press "SET" to set down and turn to next step.

4. Repeat step 2 and 3 to set the [MONTH] and [DAY].

5. Then it will enter time setting mode. There are two kinds of time format "12" hours and "24" hours. Press "MEM" button to change the time format. Press "SET" button to confirm the time format.

6. Repeat step 2 and 3 to set the [HOUR] and [MINUTE].

7. After confirming the time, the LCD will display "done" and the monitor will shut off.

Tie the Cuff

1. Plug in the cuff connecting tube into the unit.

2. Wrap the arm cuff, leaving the end of the cuff around the cuff.

3. Insert arm in the cuff. The tube must be aimed in the direction of the hand. Turn the palm of the hand upward and rest the edge of the cuff at about 2-3 cm (0.8"-1.2") over the hollow of the elbow. Tighten the cuff by pulling the end.

4. Press the hook and loop material together securely. Allow room for 2 fingers to fit between the cuff and your arm. Position the artery mark over the main artery (on the inside of your arm). **Note:** Locate the main artery by pressing with 2 fingers approximately 2 cm above the bend of your elbow on the inside of your left arm. Identify where the pulse can be felt the strongest. This is your main artery.

5. Lay your arm on a table (palm upward) so the cuff is at the same height as your heart. Make sure the tube is not kinked.

Helpful Tips

Here are a few helpful tips to help you obtain more accurate readings:

- Do not measure your blood pressure immediately after consuming a large meal.
- To obtain more accurate readings, please wait one hour after meal before measuring.
- Do not smoke or drink tea, coffee or alcohol before measuring your blood pressure.
- Do not measure your blood pressure immediately after taking a bath; wait at least 20 minutes.
- You should not be physically tired or exhausted while taking a measurement.
- Do not measure your blood pressure when you want to discharge urine.
- It is important that you relax during measurement.
- Do not take measurements if you are under stress or tense.
- Take your blood pressure at normal body temperature. If you are feeling cold or hot, wait a while before taking a measurement.
- Rest for 5 minutes before measuring.
- Wait at least 3 minutes between measurements. This allows your blood circulation to recover.
- Take the measurement in a silent room.
- Do not move and do not talk when taking a measurement.
- Do not cross your legs and keep your feet on the ground.
- Keep your back against the backrest of the chair.
- For a meaningful comparison, try to measure under similar conditions. For example, take daily measurements at approximately the same time, on the same arm, or as directed by a physician.

Start the Measurement

1. When the monitor is off, press the to turn on the monitor and start the measurement.

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