3M™ Cavilon™ Durable Barrier Cream

Product Description
3M™ Cavilon™ Durable Barrier Cream is a pH balanced, hypoallergenic, concentrated cream that provides skin protection from body fluids. The cream protects skin by providing a long lasting barrier. The cream also moisturizes and conditions skin. Use of this cream allows adhesive products to stick to the skin.

Indications for Use
Protection of at-risk skin from damage associated with incontinence and related symptoms such as anal itch. It may also be used to moisturize and condition dry skin. This product is safe to use on both intact and injured skin for example associated with incontinence.

Ingredients:
Water (Aqua), Dimethylol Adipate, Acrylate Terpolymer, Cocox Nucifer (Coco) Oil, PPG-15 Stearyl Ether, Glycerin, Isopropyl Palmitate, White Mineral Oil (Paraffinum Liquidum), Trimethylsilyloctylsilicate, Paraffin, Dimethicone, Phenoxethanol, Benzoin Acid, Dihydroxyacetone Acid, Magnesium Sulfate.

Warnings
FOR EXTERNAL USE ONLY. Avoid contact with eyes.
If condition worsens or does not improve within 7 days, consult a physician. Keep out of reach of children. In cases of ingestion, seek medical help.
Avoid use on persons with known allergy or sensitivity to any ingredient listed.
Cavilon Durable Barrier Cream is not intended to treat anal disease.

Precautions
Cavilon Durable Barrier Cream may increase the adherence of some adhesive products.
Directions for Use
Gentle cleanse skin as needed. Apply Cavilon Durable Barrier Cream sparingly to cover entire affected area.
Repeat as necessary. If the feel is oily, too much has been applied.

Storage/ Shelf Life/Disposal
For best results, store at room temperature. Do not freeze.
Avoid excessive heat and humidity.
For shelf life, refer to the expiration date on the package.
Please report a serious incident occurring in relation to the device to 3M and the local competent authority (EU) or local regulatory authority.
For further information, please contact your local 3M representative or contact us at 3M.com and select your country.

Symbol Title	Symbol	Description and Reference
Authorized Representative in European Community		Indicates the authorized representative in the European Community. ISO 15223, 5.1.2
Batch code		Indicates the manufacturer's batch code so that the batch or lot can be identified. ISO 15223, 5.1.5
Catalogue number		Indicates the manufacturer's catalogue number so that the medical device can be identified. ISO 15223, 5.1.6
Caution		Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself. Source: ISO 15223, 5.4.4
Importer		Indicates the entity importing the medical device into the EU.
CE Mark		Indicates conformity to European Union Medical Device Regulation or Directive.
Date of Manufacture		Indicates the date when the medical device was manufactured. Source: ISO 15223, 5.1.3
Do not re-use		Indicates a medical device that is intended for one use or for use on a single patient during a single procedure. Source: ISO 15223, 5.4.2
Fragrance free		Indicates medical device does not contain a fragrance additive
Green Dot		Indicates a financial contribution to national packaging recovery company per European Directive No. 94/62 and corresponding national law. Packaging Recovery Organization Europe.
Manufacturer		Indicates the medical device manufacturer as defined in EU Directives 90/385/EEC, 93/42/EEC and 98/79/EC. Source: ISO 15223, 5.1.1
Medical Device		Indicates the item is a medical device.
Natural rubber latex is not present		Indicates natural rubber or dry natural rubber latex is not present as a material of construction within the medical device or the packaging of a medical device. Source: ISO 15223, 5.4.5 and Annex B
Non Sterile		Indicates a medical device that has not been subjected to a sterilization process. Source: ISO 15223, 5.2.7
Parabens free		Indicates medical device does not contain a paraben additive.
Use-by date		Indicates the date after which the medical device is not to be used. Source: ISO 15223, 5.1.4

3392GS is only product where "DO NOT REUSE" applies.

Made in USA with globally sourced materials

3M Company
2510 Conway Ave.
St. Paul, MN 55144 USA
1-800-225-3957 (USA only)
3M.com

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