
Thermometer User Manual

FDIR-V12

Please read the guide carefully before use and keep it well.
For American please refer to "°F",
for European please refer to "°C".

1	Product introduction	3
2	Product component	3
3	Product specification	3
4	Intended use and application	4
5	Contraindication	4
6	Attention	4
7	Installation and instruction	5
8	Knowledge of body temperature	8
9	Calibration	8
10	Trouble-shooting	8
11	Cleaning instructions and product maintenance	9
12	Final disposal	9
13	Contents and Accessories	10
14	Explanation of standardized symbol	10
15	Electromagnetic compatibility information	10
16	Standard list	13
17	Right of Complaint	14

User manual version: V1.0
Issued:2023/02/16

No response	The batteries are dead? The batteries type are wrong or reversed polarity? Poor battery connection?	Replace batteries with new ones. Remove the batteries and adjust them correctly.
	Temperature is measured beyond temperature range from 32.0°C ~42.9°C (89.6°F ~109.2°F).	Refer to the user manual and measure the temperature again.
	Hi: means that the temperature measured is too high Lo: means that the temperature measured is too low	
	The ambient temperature is beyond the operating range of 50°F ~104°F (10°C ~40°C)	Keep the thermometer in the environment where the temperature is between 50°F ~104°F (10°C ~40°C for 30) minutes.
	Battery is running low	Please replace the battery in time not to disturb usage.
	The battery is almost dead and can not be used anymore.	Must replace new battery to continue usage.
	Hardware damage	Please contact distributor for support

11. Please contact distributor for support

- The probe is the most delicate part of the thermometer. It has to be clean and intact to ensure accurate readings. In order to clean the probe, please follow the way shows below: Wipe the surface of the probe gently using cotton or soft cloth with alcohol (75% Isopropyl).
- Use soft and dry cloth to clean the screen and thermometer casing. If the thermometer casing is too dirty, wipe it with soft cloth along with alcohol.
- The device is not water-proofed, please do not use detergent or place it in water or other liquid.
- We do not authorize any institution or individual to maintain and repair the product. If you suspect that the device might have issues concerning its functions, please do not repair the thermometer by yourselves.
- Thermometer is a very precise product, any improper repair or disassemble will cause inaccuracy of the measuring results.
- Please contact the after-sales service and support for any product issues within the warranty period.

10. Trouble-shooting

Problem or error message	Checklists	Solution
8		


1. Product introduction

Thanks so much for choosing our product. This high-tech infrared thermometer is used to measure the surface temperature of the forehead and convert it to the actual human temperature by an algorithm (adjustment mode). It can help you learn about the health of you and your family quickly at any time and anywhere.


Product name: Infrared thermometer
Product model: FDIR-V12

2. Product component

The thermometer mainly consists of plastic shell, infrared temperature sensor, PCBA, buzzer chip, LCD display and battery. Product front and back view shows below:

LCD display description

3. Product specification

Model	FDIR-V12
Power supply	DC 3V, 2*1.5V AAA batteries
Measuring range	Measurement Range for Body: 32.0°C~42.9°C (89.6°F~109.2°F) Measurement Range for Object Surface: 0°C~60.0°C (32.0°F~140.0°F)
Measuring accuracy	Body Mode: ±0.2°C/0.4°F, if temperature at 35.0°C ~42.0°C ±0.3°C/0.5°F, if temperature lower than 35.0°C or higher than 42.0°C Object Mode: ±1°C
Memory function	19 groups of data memory
Display resolution	0.1°C/0.1°F
Clinical repeatability	±0.3°C/0.5°F
Measuring position	Forehead
Operation condition temperature	Temperature: 50.0°F~104.0°F (10.0°C~40.0°C) Relative humidity: ≤95%RH Atmospheric pressure: 70 kPa~106 kPa

12. Final disposal

Please do not dispose of the product in the household waste at the end of its useful life. To protect the environment, dispose of empty batteries at appropriate collection sites according to national or local regulations.

13. Contents and Accessories

We only use original accessories. Check the accessory list below to ensure if the package delivered is complete.

Quantity	Parts
1pcs	FDIR-V12 device
1pcs	User Manual

14. Explanation of standardized symbol

	Complies with the European Medical Device Directive (93/42/EEC), Notified Body is SGS Belgium NV.
	Authorized representative in the European Community.
	Attention: see Instructions for use!
	Caution! Consult accompanying documents.
	Type BF applied parts
	Batch code
	Serial number
	Manufacturer information:
	IP code of the device: this device's grade of against ingress of solid foreign objects
	Disposal in accordance with Directive 2002/96/EC (WEEE)

15. Electromagnetic compatibility information

Guidance and manufacturer's declaration – electromagnetic emission – for all EQUIPMENT AND SYSTEMS

Guidance and manufacturer's declaration – electromagnetic emission
The FDIR-V12 is intended for use in the electromagnetic environment specified below. The customer or the user of FDIR-V12 should assure that it is used in such an environment.

Transport/Storage condition temperature	Temperature: -13.0°F~131.0°F (-25.0°C~55.0°C) Humidity: ≤95%RH, non-condensing Atmospheric pressure: 70 kPa~106 kPa
Grade of waterproof	IP22
Electric shock	Internally powered ME equipment
Applied part	Type BF applied part, including the whole unit
Battery life	2 years/1000 measurements
Product size	145mm*41mm*49mm
Product weight	93g
Service life	5 years
Warranty	2 years
Software version	V1.0

4. Intended use and application

This product mainly adopts infrared temperature measurement method to measure the temperature of forehead and can be used for newborns, children and adults. We recommend adults to operate the thermometer instead of babies and children.

5. Contraindication

Measuring inflammation, trauma or postoperative lesions.

6. Attention

- #### 6.1. related to measurement
- The measurement results are for reference only. Please do not make self-diagnosis and treatment according to the measurement results. If necessary, please go to hospital to take a medical treatment.
 - There is no absolute standard temperature of human body. In order to make a correct judgment for fever, it is important to know your normal body temperature, which is helpful to judge whether you have a fever or not.
 - Before measuring the forehead temperature, please make sure there is no sweat, cosmetics, oil stain, etc.
 - Before measuring, please make sure that the person does not have shower, exercise or eat in past 30 minutes. When human body is in a stable status, the body temperature measured is more referential.
 - Please do not measure the temperature near inflammation or scar, which will affect the temperature measurement results.
 - Please do not measure the body temperature immediately after taking the medicine. The temperature measured at this time is not referential.
 - Please do not measure in the environment where the temperature changes rapidly, such as the air outlet of air conditioner or heater, which will affect the temperature

Emissions test	Compliance	Electromagnetic environment- guidance
RF emissions CISPR 11	Group 1	The FDIR-V12 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The FDIR-V12 suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations flicker emissions IEC 61000-3-3	Complies	User Manual

Guidance and manufacturer's declaration – electromagnetic immunity – for all EQUIPMENT and SYSTEMS

Guidance and manufacturer's declaration – electromagnetic immunity			
The FDIR-V12 is intended for use in the electromagnetic environment specified below. The customer or the user of the FDIR-V12 should assure that it is used in such an environment.			
Immunity test	IEC 60601	Compliance level	Electromagnetic environment- guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrostatic transient / burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV differential mode ± 1 kV differential mode ± 2 kV common mode	± 1 kV differential mode ± 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% U _n ; 0.5 cycle U _n at 0°; 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% U _n ; 1 cycle and 70% U _n ; 25/30 cycles Single phase: at 0° 0% U _n ; 250/300 cycle	0% U _n ; 0.5 cycle U _n at 0°; 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% U _n ; 1 cycle and 70% U _n ; 25/30 cycles Single phase: at 0° 0% U _n ; 250/300 cycle	Mains power quality should be that of a typical commercial or hospital environment. If the user of the FDIR-V12 requires continued operation during power mains interruptions, it is recommended that the FDIR-V12 be powered from an uninterruptible power supply or a battery.

measurement results.

- When measuring repeatedly, the measurement results may have small deviation, which is a normal phenomenon.
- Please do not measure in strong electromagnetic interference environment (near working microwave oven, induction cooker, or near a active mobile phone call, etc.), which may lead to measuring improperly or inaccurate measurement results.
- This product is for personal use. In order to avoid cross infection, please pay attention to the cleaning and disinfection of the product.

6.2. related to product

- This product belongs to precision equipment. Please put it in the packaging box after use to avoid liquid splashing into the device and probe and prevent small foreign matter (such as dust) from falling into the probe, which may affect measurement results.
- Avoid falling on the ground or impacted by external force. Please do not disassemble and assemble by yourself.
- Avoid touching the probe directly with your finger or blowing it with your mouth. When the infrared probe is damaged or dirty, the measurement results may be inaccurate.
- Please put this product out of reach of the child to prevent the child from swallowing the battery or small parts of the thermometer.
- Do not put the thermometer and battery into the fire to prevent explosion.
- Please take out the batteries if the thermometer will not be used for more than 3 months.
- Self-diagnosis and treatment according to the measurement results are dangerous. Please consult a professional doctor for treatment based on the measurement results.

7. Installation and instruction


7.1. Installation of product

Put two AAA-batteries into the battery house on the back of the device. At the point the product will start self-inspection, and then turn into power off automatically. (If the battery power is low when starting up, please replace the battery to ensure adequate power supply)

7.2. Measurement process

7.2.1. Forehead temperature measurement

- Put the probe toward the forehead directly at a distance from 1cm to 6cm.


Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Note: U _i is the a.c. mains voltage prior to application of the test level.			

Guidance and manufacturer's declaration – electromagnetic immunity – for EQUIPMENT and SYSTEM

Guidance and manufacturer's declaration – electromagnetic immunity			
The FDIR-V12 is intended for use in the electromagnetic environment specified below. The customer or the user of the FDIR-V12 should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment- guidance
3 Wms		3V	Portable and mobile RF communications equipment should be used no closer to any part of the FDIR-V12, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = \left[\frac{3.5}{U_i} \right] \sqrt{P}$ $d = \left[\frac{12}{U_i} \right] \sqrt{P}$
150 kHz to 80 MHz		150 kHz to 80 MHz	
6 V in ISM and amateur radio bands between 0.15 MHz and 80 MHz		6 V in ISM and amateur radio bands between 0.15 MHz and 80 MHz	
Conducted RF IEC 61000-4-6	10 V/m	10 V/m	
80 MHz to 2.7 GHz		80 MHz to 2.7 GHz	
Radiated RF IEC 61000-4-3	385MHz; 5785MHz Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communication equipment (Refer to table 9 of IEC 60601-1-2:2014)	385MHz; 5785MHz Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communication equipment (Refer to table 9 of IEC 60601-1-2:2014)	Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol:

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.
Note 2: These guidelines may not apply in all situations.
Electromagnetic is affected by absorption and reflection from structures, objects and people.

- Press the [Power/Scan Button] lightly. Testing result can be shown on display within 1 second.


7.2.2. Object Temperature Measurement

- Press the [Mode Switch Button] slightly for object mode.
- Put the probe towards the object.
a. At distance 1-5 cm in a perpendicular a line.
b. When measuring the temperature of food or drink, be careful not to measure on the packaging but close to the contents themselves.

- Press the [Power/Scan Button] lightly and testing result can be shown on display within 1 second.

7.2.3. Sound and backlight explanation

If infrared thermometer beeps once, it means everything goes well. If the measurement result is 37.6°C or above, the infrared thermometer will beep four times. Measurement readings will be displayed on screen.

Measurement range for body	Buzzer warning	Color of backlight
32.0°C (89.6°F) ≤ T ≤ 37.5°C (99.6°F)	The buzzer warning by sending out a long "beep"	Green
37.6°C (99.7°F) ≤ T ≤ 42.9°C (109.2°F)	The buzzer warning by sending out a "beep" "beep" "beep" "beep" "beep" "beep" "beep" "beep"	Red
T < 32.0°C (89.6°F) Displaying Lo	The buzzer warning by sending out a "beep" "beep" "beep" "beep"	Green
T > 42.9°C (109.2°F) Displaying Hi	The buzzer warning by sending out a "beep" "beep" "beep" "beep"	Green

Measurement range for object	Buzzer warning	Color of backlight
0°C(32.0°F) ≤ T ≤ 60.0°C (140.0°F)	The buzzer warning by sending out a long "beep"	Green
T < 0°C(32.0°F) Displaying Lo	The buzzer warning by sending out a "beep" "beep" "beep"	Green
T > 60.0°C (140.0°F) Displaying Hi	The buzzer warning by sending out a "beep" "beep" "beep"	Green

- The ISM (industrial, scientific and medical) bands between 150 kHz and 80 MHz are 6,765 MHz to 6,795 MHz, 13,553 MHz to 13,567 MHz, 26,957 MHz to 27,283 MHz and 40,66 MHz to 40,70 MHz. The amateur radio bands between 0,15 MHz and 80 MHz are 1,8 MHz to 2,0 MHz, 3,5 MHz to 4,0 MHz, 5,3 MHz to 5,4 MHz, 7 MHz to 7,3 MHz, 10,1 MHz to 10,15 MHz, 14 MHz to 14,2 MHz, 18,07 MHz to 18,17 MHz, 21,0 MHz to 21,4 MHz, 24,89 MHz to 24,99 MHz, 28,0 MHz to 29,7 MHz and 50,0 MHz to 54,0 MHz.
- Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless), telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the FDIR-V12 is used exceeds the applicable RF compliance level above, the FDIR-V12 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the FDIR-V12.
- Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3V/m.

Recommended separation distances between portable and mobile RF communications equipment and the EQUIPMENT or SYSTEM - for EQUIPMENT and SYSTEMS

Recommended separation distances between portable and mobile RF communications equipment and the FDIR-V12			
The FDIR-V12 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the FDIR-V12 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the FDIR-V12 as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz outside ISM and amateur radio bands $d = \left[\frac{12}{U_i} \right] \sqrt{P}$	150 kHz to 80 MHz in ISM and amateur radio bands $d = \left[\frac{3.5}{U_i} \right] \sqrt{P}$	800 MHz to 2.7 GHz $d = \left[\frac{1}{U_i} \right] \sqrt{P}$
0.01	0.12	0.20	0.035
0.1	0.38	0.63	0.11
1	1.2	2.00	0.35
10	3.8	6.32	1.10
100	12	20.00	3.5
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol:			

16. Standard list

FDIR-V12 complies with the following standards

EN ISO 15223-1	Medical device – symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements
EN 1041	Information supplied by the manufacturer with medical devices
IEC 60601-1	Medical electrical equipment Part 1: General requirements for basic safety and essential performance

7.2.4. Memory function

This product has 19 groups of data memory function. After power on, press the [Memory/Voice Button] shortly to query the user's historical test data. During display process, the current data group sequence is displayed first, and then comes the temperature data of this group sequence.

2. Continue to press [Memory/Voice Button] shortly to query the next group of data. If it exceeds the last group, it will display the first group data again.

3. If there is no operation for 3 seconds, it will exit the memory mode display.

7.2.5. Mute mode function

After power on, press [Memory/Voice Button] for about 2 seconds to turn on or off the mute mode.

7.2.6. Switching measurement modes
Once power on, object mode and body mode can be switched by pressing [Mode Switch Button] slightly.

7.2.7. Power off
If the device is not in use, it will Power off automatically in 30 seconds.

7.3. Temperature unit selection

A. After power on, press and hold the [Mode Switch Button] for 10 seconds. The measurement system (°C or °F) automatically appears on the display.

B. Then still keep holding the [Mode Switch Button] for 2 more seconds, the unit °C and °F will be switched automatically. Release the button when the selected unit appeared. The device is auto into measure ready status. The unit selected is active.

IEC 60601-1-11	Medical electrical equipment – Part 1-11: General requirements for basic safety and essential performance – Collateral Standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment
IEC 60601-1-2	Medical electrical equipment – Part 1-2: General requirements for basic safety and essential performance – Collateral standard: Electromagnetic compatibility – Requirements and tests
ASTM E 1965	Standard Specification for Infrared Thermometer for Intermittent Determination of Patient Temperature
IEC 62304	Medical device software – Software life-cycle processes
IEC 62366	Medical devices – Application of usability engineering to medical devices IEC 62366:2007
ISO 80601-2-56	Medical electrical equipment – Part 2-56: Particular requirements for basic safety and essential performance of clinical thermometers for body temperature measurement
ISO 10993-1	Biological evaluation of medical devices – Part 1: Evaluation and testing within a risk management process

17. Right of Complaint

We provide 2 years of Right of Complaint starting from the date of purchase. Please refer to the followings situations that are excluded from the free repair services within the warranty period.

- All damages caused by disassembly and repair of the device by yourselves.
 - All damages caused by dropping the device during usage, or transport.
 - All damages caused by improper usage of the device and not following the instructions on the user manual.
- Please contact after-sales service and support and enclose your product purchase receipt while claiming for warranty services.

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