EN
This medical device and Instructions for Use are intended for use by a lay person. Read the instruction for use before starting the treatment with Sorbact® Absorption. Device Description Sorbact® Absorption is a bacteria and fungi binding wound dressing, based on Sorbact® Technology. It consists of a Sorbact® DACC™-coated wound contact layer combined with an absorbent core. The product absorbs and retains exudate. Sorbact® Absorption prevents and treats wound infections and facilitates the wound healing process. Sorbact® wound contact layer is shown to bind wound related pathogens up to seven days <i>in vitro</i>. Intended Purpose Sorbact® Absorption is intended for use on superficial exuding wounds such as surgical wounds and acute wounds. Who can use Sorbact® Absorption? Sorbact® Absorption is intended to be used on children and adults. The dressing can be used by healthcare professionals and lay persons in healthcare facilities and home environment. How does Sorbact® Absorption work? Sorbact® Absorption binds common microorganisms, such as <i>Staphylococcus aureus</i> (including MRSA), <i>Streptococcus</i> spp., <i>Escherichia coli</i>, <i>Pseudomonas aeruginosa</i> and <i>Candida albicans</i>, as shown <i>in vitro</i>. There is no release of antimicrobial agents to the wound. When shall Sorbact® Absorption not be used? There are no known contraindications to the use of Sorbact® Absorption. Warnings and precautions should be noted. Warnings and Precautions Do not use if you have previously experienced hypersensitivity to the dressing (dialkylcarbamoyl chloride [DACC™], cellulose acetat, cellulose, cotton, polyester and viscose). Do not use in combination with fatty products, such as ointments, creams and solutions as they may decrease the binding of bacteria. Do not re-use. Sorbact® Absorption is for single patient and single use only. Re-use may lead to cross contamination and increased risk of infection. Do not use if the peel pouch is already opened or damaged, since the sterility cannot then be guaranteed. Consult a healthcare professional if: • Wound fails to improve or deteriorates • Side effect is observed • Signs of infection such as red, warm, swollen or painful skin appear How to use Sorbact® Absorption 1. Clean the wound and dry the surrounding skin. 2. Open the packaging without touching the green side of the dressing. 3. Apply the dressing with the green side of the dressing in direct contact with the wound surface. The dressing should overlap the wound margins by at least 2 cm. Do not cut the dressing. 4. Fixate with an elastic bandage or skin friendly tape. 5. The dressing change frequency depends on the condition of the wound and the exudate level. Change the dressing before it is fully saturated or at signs of leakage. The dressing should not be left in place for more than a few days. 6. Should the dressing adhere to the wound, moisten the dressing with water or saline solution to assist removal and to avoid disruption of the healing wound. Storage and Disposal Sorbact® Absorption shall be stored dry and kept away from sunlight. Disposal should be made according to local environmental procedures. Notice to User Any serious incident that has occurred in relation to the device should be reported to ABIGO Medical AB (complaints@abigo.com) and the competent authority of your state.

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