Chlamydia Rapid Test
Package Insert
For self-testing

REF	ICH-N502H-01/ICH-N502H-02/ICH-N502H-05	English
	ICH-N502FH-01/ICH-N502FH-02/ICH-N502FH-05	
	ICH-N502MH-01/ICH-N502MH-02/ICH-N502MH-05	

【INTENDED USE】

The test is a rapid chromatographic immunoassay for the qualitative detection of *Chlamydia trachomatis* antigens in female cervical swab or male urethral swab specimens to aid in the diagnosis of Chlamydia infection.

The test is for self-testing use. The test is not automated. For *in vitro* diagnostic use only.

【SUMMARY】

Chlamydia trachomatis is the most common cause of sexually transmitted venereal infection in the world. It is composed of elementary bodies (the infectious form) and reticulate or inclusion bodies (the replicating form). *Chlamydia trachomatis* has both a high prevalence and asymptomatic carriage rate, with frequent serious complications in both women and neonates. Complications of Chlamydia infection in women include cervicitis, urethritis, endometritis, pelvic inflammatory disease (PID) and increased incidence of ectopic pregnancy and infertility.¹ Vertical transmission of the disease during parturition from mother to neonate can result in inclusion conjunctivitis or pneumonia. In men, complication of Chlamydia includes urethritis and epididymitis. At least 40% of the nongonococcal urethritis cases are associated with Chlamydia infection. Approximately 70% of women with end cervical infections and up to 50% of men with urethral infections are asymptomatic. Traditionally, Chlamydia infection has been diagnosed by detection of Chlamydia inclusions in tissue culture cells. Culture method is the most sensitive and specific laboratory method, but it is labor intensive, expensive, long (18-72 hours) and not routinely available in most situations.^{2,3}

【PRINCIPLE】

The Chlamydia Rapid Test is a lateral flow immunoassay based on the principle of the double-antibody sandwich technique. The test contains an antibody to *Chlamydia trachomatis* that is conjugated with latex particles on the conjugate pad and antibody specific to the *Chlamydia trachomatis* antigens on the test line region. During testing, the extracted antigen solution reacts with an antibody to *Chlamydia trachomatis* conjugated with latex particles to form an antibody-antigen complex. The complex migrates up to react with the antibody to *Chlamydia trachomatis* on the membrane and generates a colored line in the test region. The presence of this colored line in the test line region indicates a positive result, while its absence indicates a negative result. To serve as a procedural control, a colored line will always appear in the control line region, indicating that the test is working properly.

【REAGENTS】

The test contains *Chlamydia trachomatis* antibody conjugated latex particles on the conjugate pad and *Chlamydia trachomatis* antibody coated on the membrane.

【WARNINGS AND PRECAUTIONS】

Please read all the information in this package insert before performing the test.

- For self-testing *in vitro* diagnostic use only.
- Do not use the test after expiration date. Do not reuse the test.
- Store in a dry place at 2-30°C (36-86°F). Do not use it if any of the test kit content or packaging is damaged or opened.
- Do not eat, drink or smoke in the area where the specimens or kits are handled.
- Wash hands thoroughly before and after handling.
- Handle all specimens as if they contain infectious agents. The used tests and specimens may be infectious and should be discarded according to local regulations.
- Use the test only once. Do not dismantle and touch the test window of the test cassette.
- The kit must not be frozen or used after the expiration date printed on the package.
- Keep out of the reach of children.
- This product is not intended for use with children.
- Carefully handle the buffer - do not ingest and avoid contact with skin or eyes. If contact does occur, immediately rinse with running water.
- Do not swallow Buffer 1, Buffer 2 or any other components in the box.
- The buffer 1 contains 0.8% Sodium hydroxide. IF ON SKIN: Wash with plenty of soap and water. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If skin or eye irritation persists: Get medical advice/attention.
- It is important that the timings are strictly followed to ensure accurate results.
- Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority.
- Components provided in the kit are approved for use in the Chlamydia Rapid Test. Do not use any other commercial kit component.
- The test is not a substitute for professional medical advice. Do not take any decision of medical relevance without first consulting the appropriate healthcare professional.
- The swab may be initially uncomfortable, but this is perfectly normal.

Warning

【STORAGE AND STABILITY】

The kit can be stored at room temperature or refrigerated (2-30°C). The test is stable through the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. **DO NOT FREEZE.** Do not use after the expiration date.

Note: It is suggested to use test cassette within one hour after removing it from the foil pouch.

【MATERIALS】

REF	Components						Tube holder
	Test cassette	Package insert	Sterile female cervical swab	Sterile male urethral swab	Buffer 1	Buffer 2	
ICH-N502H-01	1	1	1	1	1	1	On the box
ICH-N502H-02	2	1	2	2	2	2	
ICH-N502H-05	5	1	5	5	5	5	
ICH-N502FH-01	1	1	1	/	1	1	
ICH-N502FH-02	2	1	2	/	2	2	
ICH-N502FH-05	5	1	5	/	5	5	
ICH-N502MH-01	1	1	/	1	1	1	
ICH-N502MH-02	2	1	/	2	2	2	
ICH-N502MH-05	5	1	/	5	5	5	

*Note:

Buffer 1 0.3 mL (0.8% NaOH)		Buffer 2: 0.28 mL (0.73% HCl), yellow
---------------------------------------	--	---

Materials Required But Not Provided

- Timer

【COLLECTING THE SPECIMEN AND DIRECTIONS FOR USE】


Before performing the test, swab specimen must be collected following the instructions below.

1. Wash your hands with soap and rinse with clear water.

To collect **Male Urethral Swab Specimens:**


Do not urinate for at least 1 hour prior to specimen collection.

Remove the sterile male urethral swab from the pouch, insert the male urethral swab tip into the urethra about 2-4 cm, rotate the swab 360° in one direction (clockwise or counter clockwise), let stand for 10 seconds, and then withdraw. Carefully remove the swab avoiding contact with any external surfaces.

To collect **Female Cervical Swab Specimens:**

Remove the sterile female cervical swabs from the pouch. Insert the swab tip into vaginal tract for half a minute and then carefully remove the swab avoiding contact with any external surfaces.


Note: Specimen should be tested immediately.


- Remove the cover of the Buffer 1 tube and insert the tube into the tube holder on the packaging box. Buffer 1 is colorless. Immediately insert the swab into the tube with Buffer 1, compress the bottom of tube and rotate swab 15 times. Let stand for 2 minutes.
- Break off the Buffer 2 at the breaking point (the narrowest place). Then add entire Buffer 2 solution to the Buffer 1 tube. The solution would turn turbid. Compress the bottle of tube and rotate the swab 15 times until the solution turn clear with a slight green or blue tint. If the swab is bloody, the color will turn yellow or brown. Let stand 1

minute.

- Press the swab against the side of tube and withdraw the swab while squeezing the tube. Keep as much liquid in the tube as possible. Close the cap onto the tube, then open the small cap.

- Remove the test cassette from the sealed foil pouch and use it within one hour. Place the test cassette on a clean and level surface.
- Invert the Buffer tube and add **3 full drops of the extracted solution** to the specimen well(S) of the test cassette, and then start the timer. Avoid trapping air bubbles in the specimen well. Do not move the test cassette during test developing.
- Wait for the color to appear. **Read the result at 10 minutes;** do not interpret the result after 20 minutes.

【READING THE RESULTS】

POSITIVE: Two colored lines appear. Both T (Test) line and C (Control) line appear.

A positive result indicates that *Chlamydia trachomatis* antigens were detected in the specimen.

*NOTE: The intensity of the color in the test line region (T) will vary depending on the concentration of Chlamydia present in the specimen.

NEGATIVE: One colored line appears in the control line region (C). No line appears in the test line region (T).

A negative result indicates that *Chlamydia trachomatis* antigens are not present in the specimen, or are present below the detectable level of the test.

INVALID: Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

【QUALITY CONTROL】

A procedural control is included in the test. A colored line appearing in the control line region(C) is considered an internal procedural control. It confirms the test is working properly.

【PERFORMANCE CHARACTERISTICS】

Sensitivity	Specificity
96.0%	98.3%

Performance data obtained by clinical study with 220 participants enrolled. PCR has been utilized as reference method.

Sensitivity

The Chlamydia Rapid Test can detect the level as low as 5E+06 CFU/mL(5E+04CFU/test) for Chlamydia.

Hook Effect

No high-dose hook effect was detected at the concentration up to 2.5E+09 CFU/mL.

Precision-Repeatability

No difference was detected in intra lot.

Precision-Reproducibility

No difference was detected between days, sites, lots and operators.

Cross-Reactivity

No cross-reactivity was observed when tested with following organisms:
Acinetobacter calcoaceticus *Pseudomonas aeruginosa*
Proteus mirabilis *Acinetobacter spp*
Neisseria meningitidis *Neisseria gonorrhoeae*
Enterococcus faecalis *Salmonella choleraesuis*
Group B/C Streptococcus *Candida albicans*
Hemophilus influenzae *Staphylococcus aureus*
Proteus vulgaris *Branhamella catarrhalis*
Klebsiella pneumoniae *Gardnerella vaginalis*
Chlamydia Pneumoniae

Interfering Substances

The following potentially interfering substances do not interfere with the test results.

Name	Name
Nifuratel and Nystatin Vaginal suppositories	Compound Metronidazole Vaginal Suppositories
Amphotericin B Vaginal Effervescent Tablets	Femfresh Daily Intimate Wash
Terbinafine Hydrochloride Vaginal	eve™ Feminine Wash
Nifuratel Tablets	Miconazole Nitrate Vaginal Suppositories
Azithromycin	Ofloxacin
Erythromycin	Amoxicillin

【LIMITATIONS】

1. The Chlamydia Rapid Test is for *in vitro* diagnostic use only. This test should be used for the detection of *Chlamydia trachomatis* antigens from female cervical swab or male urethral swab specimens. Neither the quantitative value nor the rate of increase in *Chlamydia trachomatis* antigens concentration can be determined by this qualitative test.
2. This test will only indicate the presence of *Chlamydia trachomatis* antigens in specimens from both viable and non-viable Chlamydia. Performance with specimens other than female cervical swabs and male urethral swabs has not been assessed.
3. Therapeutic failure or success cannot be determined as antigen may persist following appropriate antimicrobial therapy.
4. Excessive blood on the swab may cause false positive results.
5. Endocervical specimens from females should not be collected during menstrual period.

【EXTRA INFORMATION】

1. How does the Chlamydia test cassette work?

Chlamydia trachomatis is the most common cause of sexually transmitted venereal infection in the world. The Chlamydia Rapid Test detects specifically the antigens in vaginal discharge or urethral discharge to ascertain the presence of the chlamydia.

2. When should the test be used?

Chlamydia trachomatis antigens can be detected from genital infections, and it is recommended to conduct tests in individuals with symptoms who meet the case definition (vaginitis, urethritis, epididymitis, prostatitis, or the following symptoms: dysuria, urethral discomfort, pain, swelling, and tenderness of the epididymis, erythema, redness, and edema of the scrotal skin, swelling, erythema, severe pain, and thickening of the vas deferens, lumbago, lower abdominal heaviness, bleeding, diarrhea, and mucoid discharge, abnormal vaginal discharge, bleeding outside of menstruation or after intercourse, and symptoms such as dysuria, frequent urination, and urgency).and to test asymptomatic individuals limited to contacts of confirmed Chlamydia cases or probable cases and to at-risk health workers.

3. Can the result be incorrect?

The results are accurate as far as the instructions are carefully respected. Nevertheless, the result can be incorrect if Chlamydia Rapid Test gets wet before performing the test or if the quantity of specimens dispensed in the specimen well is too much or not sufficient, or if the number of extracted specimens drops are less than 3 or more than 4.

Detection of Chlamydia is dependent on the number of organisms present in the specimen. This can be affected by specimen collection methods and patient factors such as age, history of Sexually Transmitted Diseases (STDs), presence of symptoms, etc. The minimum detection level of this test may vary according to serovar.

Besides, due to immunological principles involved, there exist the chances of false results in rare cases. A consultation with the doctor is always recommended for such tests based on immunological principles.

4. How to interpret the test if the color and the intensity of the lines are different?

The colour and intensity of the lines have no importance for result interpretation. The lines should only be homogeneous and clearly visible. The test should be considered as positive whatever the colour intensity of the test line is.

5. What is the line that appears under the in the control line region (C) for?

When this line appears, it only means that the test unit is performing well.

6. What do I have to do if the result is positive?

If the result is positive, it means that the Chlamydia was detected in vaginal discharge or urethral discharge and that you should consult a doctor to show the test result. Then, the doctor will decide whether additional analysis should be performed.




7. What do I have to do if the result is negative?

If the result is negative, it means that it was not possible to detect the *Chlamydia trachomatis* antigens. However, if the symptoms persist, it is recommended to consult a physician.

【BIBLIOGRAPHY】

1. Sanders J.W. et al. Evaluation of an Enzyme Immunoassay for Detection of Chlamydia trachomatis in Urine of Asymptomatic Men. J.Clinical Microbiology, (1994) 32, 24-27.
2. Jaschek, G. et al. Direct Detection of *Chlamydia trachomatis* in Urine Specimens from Symptomatic and Asymptomatic Men by Using a Rapid Polymerase Chain Reaction Assay. J. Clinical Microbiology, (1993) 31, 1209-1212.
3. Schachter, J. Sexually transmitted *Chlamydia trachomatis* infection. Postgraduate Medicine, (1982) 72, 60-69.

Index of Symbols

	Consult instructions for use or consult electronic instructions for use		Contains sufficient for <n> tests		Temperature limit
	<i>In vitro</i> diagnostic medical device		Batch code		Catalogue number
	Authorized representative in the European Community/European Union		Use-by date		Do not re-use
	Do not use if package is damaged and consult instructions for use		Manufacturer		Unique device identifier
	Warning		Importer		

Hangzhou AllTest Biotech Co.,Ltd.

#550,Yinhai Street
Hangzhou Economic & Technological Development Area
Hangzhou, 310018 P.R. China
Web: www.alltests.com.cn Email: info@alltests.com.cn

EU REP

MedNet EC-REP GmbH
Borkstrasse 10,
48163 Muenster,
Germany

Summary of safety and performance is available on website:
<https://ec.europa.eu/tools/eudamed/>

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【BIBLIOGRAFI】

1. Sanders J.W. et al. Evaluation of an Enzyme Immunoassay for Detection of *Chlamydia trachomatis* in Urine of Asymptomatic Men. J.Clinical Microbiology, (1994) 32, 24-27.
2. Jaschek, G. et al. Direct Detection of *Chlamydia trachomatis* in Urine Specimens from Symptomatic and Asymptomatic Men by Using a Rapid Polymerase Chain Reaction Assay. J. Clinical Microbiology, (1993) 31, 1209-1212.
3. Schachter, J. Sexually transmitted *Chlamydia trachomatis* infection. Postgraduate Medicine, (1982) 72, 60-69.

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