Front Panel
When Folded

Back Panel
When Folded

- de Pen-Nadeln
- es Aghi penna
- es Agujas para pluma de inyección
- es Jehly pro injekční pera
- es Βελόνες για συσκευή τύπου πέννας
- sv Pennkanyler
- fi Kynäneulat
- de Penkanyler

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NOTES:

- Dimensions: Flat: Width: 220mm. Height: 440mm
Folded: Width: 55mm. Height: 110mm.
- Material is referenced on MBS-PN-013 (Suzhou Manufacturing). No change may be made from the material specified without the prior approval of embecta.
- Dimensional Tolerance: +/- 1mm.
- Red lines () do not print.

 v13		Drawing No. PCI89105		Colours	
embecta Medical II LLC 300 Kimball Drive Parsippany, NJ 07054 USA		Sheet 3 of 4	SAP Rev 01	1. Process Cyan 2. Process Magenta 3. Process Yellow 4. Process Black	
Title embecta™ Pen Needles•3bevel•TW Insert					
Version Letter A		Cycle Number 06		Part No. N/A ECO No. N/A	
Date 01APR2025		Page Size 220 x 440mm			
Notes Template-32panel-Insert-V1					
The subject of this drawing is submitted in confidence and it and the drawing will and shall remain the property of embecta and are not to be used without its written permission and approval. ©2025 embecta, Parsippany, New Jersey, USA.					
For Physical Requirements and Printing Requirements refer to the NOTES section of this specification.					

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
 
2797 PAP

embecta and the embecta logo are trademarks of Embecta Corp. and registered or used in the United States or other countries. © 2025 Embecta Corp. All rights reserved.

embecta.com
20250415 00000001

The diagram illustrates the process of folding a square piece of paper into a cube. It is divided into three sections, each showing a different stage of the process:

- Fold 1:** A square piece of paper is shown. The left half is labeled "Front" and the right half is labeled "Back".
- Fold 2:** The square is folded in half, with the "Front" and "Back" labels visible on the resulting rectangle.
- Fold 3:** The folded square is folded again, with the "Front" label visible on the resulting smaller rectangle.

NOTES:

1. Dimensions: Flat: Width: 220mm. Height: 440mm
Folded: Width: 55mm. Height: 110mm.
2. Material is referenced on MBS-PN-013 (Suzhou Manufacturing).
No change may be made from the material specified without the prior approval of embecta.
3. Dimensional Tolerance: +/- 1mm.
4. Red lines () do not print.

 embecta Medical II LLC Kimball Drive Parsippany, NJ 07054 USA	Drawing No. PCI89105		Colours 1. Process Cyan 2. Process Magenta 3. Process Yellow 4. Process Black
	Sheet 4 of 4	SAP Rev 01	
embecta™ Pen Needles•3bevel•TW Insert			Part No. N/A ECO No. N/A
on Letter A	Cycle Number 06	Notes Template-32panel-Insert-V1	
01APR2025	Page Size 220 x 440mm	For Physical Requirements and Printing Requirements refer to the NOTES section of this specification.	
Subject of this drawing is submitted in confidence and it and its drawing will and shall remain the property of embecta and are not to be used without its written permission and approval. © 2025 embecta. Parsippany, New Jersey, USA.			