

OptiChamber Diamond

Anti-static valved holding chamber

Single patient use

[illegible][illegible]

Version	ECR Number	Drawn		Checked		Approved		1109825
		Initials	Date	Initials	Date	Initials	Date	
00	9835	AM	30JUL13	KB	31JUL13	TS	31JUL13	
01	10999	AM	21JUL15	DP	22 JUL 15	AA	22 JUL 15	
02	11025	KAB	07 SEP 15	<i>Julius</i>	<i>07-SEPT-15</i>	<i>/A</i>	<i>14-SEP-15</i>	

SPECIFICATIONS

Format:	Mini-fold instructions for use
Flat (trim) size:	Metric: 625.4mm (w) x 750mm (h) Imperial: 24.6" (w) x 29.52" (h)
Flat tolerance:	Metric: ± 2 mm Imperial: ± 0.078 "
Final finished size:	Metric: 78mm (w) x 125mm (h) Imperial: 3.07"(w) x 4.92"(h)
Finished tolerance:	Metric: ± 2 mm Imperial: ± 0.078 "
Finish/folding:	Trim and accordion fold to final finished size.
Stock / Coating:	30# White
Ink:	1color / 2sides. Black only.
Bar code:	GS1-128
Proof:	First article inspection is required for approval Production commencement only after approval.
Artwork:	Do not use copy for production. For reference only.

DO NOT SCALE CONFIDENTIAL:- THIS DRAWING CONTAINS PROPRIETARY INFORMATION AND SHALL NOT BE USED. REPRODUCED OR DISTRIBUTED WITHOUT THE PRIOR WRITTEN CONSENT OF RESPIRONICS	PHILIPS  RESPIRONICS Respironics Respiratory Drug Delivery (UK) Ltd Chichester Business Park, City Fields Way, Tangmere, Chichester, PO20 2FT, UK	DESCRIPTION	DOC TYPE:-
		P2 IFU DIAMOND ROW MATERIAL NUMBER 1109825	LBL Page 1 of 3